

LSP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 37kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human LSP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

WP34; pp52;

Catalog #: A5617

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4046

Isotype: IgG

Swiss Prot: P33241

Purity: Affinity purification

For research use only.

Background:

This gene encodes an intracellular F-actin binding protein. The protein is expressed in lymphocytes, neutrophils, macrophages, and endothelium and may regulate neutrophil motility, adhesion to fibrinogen matrix proteins, and transendothelial migration. Alternative splicing results in multiple transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).