

PGRMC1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 21kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PGRMC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MPR; HPR6.6;

Catalog #: A5619

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10857

Isotype: IgG

Swiss Prot: O00264

Purity: Affinity purification

For research use only.

Background:

This gene encodes a putative membrane-associated progesterone steroid receptor. The protein is expressed predominantly in the liver and kidney.

To place an order, please [Click HERE](#).