

C9

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 63kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human C9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

C9D; ARMD15;

Catalog #: A5622

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 735

Isotype: IgG

Swiss Prot: P02748

Purity: Affinity purification

For research use only.

Background:

This gene encodes the final component of the complement system. It participates in the formation of the Membrane Attack Complex (MAC). The MAC assembles on bacterial membranes to form a pore, permitting disruption of bacterial membrane organization. Mutations in this gene cause component C9 deficiency.

To place an order, please [Click HERE](#).