

PLTP

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 54kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PLTP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BPIFE; HDLCQ9;

Catalog #: A5628

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5360

Isotype: IgG

Swiss Prot: P55058

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is one of at least two lipid transfer proteins found in human plasma. The encoded protein transfers phospholipids from triglyceride-rich lipoproteins to high density lipoprotein (HDL). In addition to regulating the size of HDL particles, this protein may be involved in cholesterol metabolism. At least two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).