

## ECE1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**87kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human ECE1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ECE;

**Catalog #:**A5638

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**1889

**Isotype:**IgG

**Swiss Prot:**P42892

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is involved in proteolytic processing of endothelin precursors to biologically active peptides. Mutations in this gene are associated with Hirschsprung disease, cardiac defects and autonomic dysfunction. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

*To place an order, please [Click HERE](#).*