

IFNL3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 22kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant Protein of human IFNL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IL28B; IL28C; IL-28B;

Catalog #: A5648

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 282617

Isotype: IgG

Swiss Prot: Q8IZI9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a cytokine distantly related to type I interferons and the IL-10 family. This gene, interleukin 28A (IL28A), and interleukin 29 (IL29) are three closely related cytokine genes that form a cytokine gene cluster on a chromosomal region mapped to 19q13. Expression of the cytokines encoded by the three genes can be induced by viral infection. All three cytokines have been shown to interact with a heterodimeric class II cytokine receptor that consists of interleukin 10 receptor, beta (IL10RB) and interleukin 28 receptor, alpha (IL28RA).

To place an order, please [Click HERE](#).