

OCA2

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 92kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human OCA2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P; BEY; PED; BEY1; BEY2; BOCA; EYCL; HCL3; EYCL2; EYCL3; SHEP1; D15S12;

Background:

This gene encodes the human homologue of the mouse p (pink-eyed dilution) gene. The encoded protein is believed to be an integral membrane protein involved in small molecule transport, specifically tyrosine - a precursor of melanin. Mutations in this gene result in type 2 oculocutaneous albinism.

To place an order, please [Click HERE](#).

Catalog #: A5653

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4948

Isotype: IgG

Swiss Prot: Q04671

Purity: Affinity purification

For research use only.