

TNFRSF13B

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 32kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TNFRSF13B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CVID; RYZN; TACI; CD267; CVID2; IGAD2; TNFRSF14B;

Catalog #: A5658

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23495

Isotype: IgG

Swiss Prot: O14836

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a lymphocyte-specific member of the tumor necrosis factor (TNF) receptor superfamily. It interacts with calcium-modulator and cyclophilin ligand (CAML). The protein induces activation of the transcription factors NFAT, AP1, and NF-kappa-B and plays a crucial role in humoral immunity by interacting with a TNF ligand. This gene is located within the Smith-Magenis syndrome region on chromosome 17.

To place an order, please [Click HERE](#).