

## MME

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IP FC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IP 1:20 - 1:50 FC 1:20 - 1:50

**Calculated MW:** 85kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human MME

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

i

**Synonym:**

NEP; SFE; CD10; CALLA;

**Catalog #:** A5664

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 4311

**Isotype:** IgG

**Swiss Prot:** P08473

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a common acute lymphocytic leukemia antigen that is an important cell surface marker in the diagnosis of human acute lymphocytic leukemia (ALL). This protein is present on leukemic cells of pre-B phenotype, which represent 85% of cases of ALL. This protein is not restricted to leukemic cells, however, and is found on a variety of normal tissues. It is a glycoprotein that is particularly abundant in kidney, where it is present on the brush border of proximal tubules and on glomerular epithelium. The protein is a neutral endopeptidase that cleaves peptides at the amino side of hydrophobic residues and inactivates several peptide hormones including glucagon, enkephalins, substance P, neurotensin, oxytocin, and bradykinin. This gene, which encodes a 100-kD type II transmembrane glycoprotein, exists in a single copy of greater than 45 kb. The 5' untranslated region of this gene is alternatively spliced, resulting in four separate mRNA transcripts. The coding region is not affected by alternative splicing.

**To place an order, please [Click HERE](#).**