

IL3

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 17kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IL-3; MCGF; MULTI-CSF;

Catalog #: A5671

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3562

Isotype: IgG

Swiss Prot: P08700

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a potent growth promoting cytokine. This cytokine is capable of supporting the proliferation of a broad range of hematopoietic cell types. It is involved in a variety of cell activities such as cell growth, differentiation and apoptosis. This cytokine has been shown to also possess neurotrophic activity, and it may be associated with neurologic disorders.

To place an order, please [Click HERE](#).