

CD58

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 28kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CD58

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ag3; LFA3; LFA-3;

Catalog #: A5684

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 965

Isotype: IgG

Swiss Prot: P19256

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the immunoglobulin superfamily. The encoded protein is a ligand of the T lymphocyte CD2 protein, and functions in adhesion and activation of T lymphocytes. The protein is localized to the plasma membrane. Alternatively spliced transcript variants have been described.

To place an order, please [Click HERE](#).