

SH3GL2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50- 1:200 IF 1:50 - 1:200

Calculated MW: 39kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SH3GL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CNSA2; SH3P4; EEN-B1; SH3D2A; Endophilin-A1

Catalog #: A5701

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6456

Isotype: IgG

Swiss Prot: Q99962

Purity: Affinity purification

For research use only.

Background:

Implicated in synaptic vesicle endocytosis. May recruit other proteins to membranes with high curvature.

To place an order, please [Click HERE](#).