

FOXP3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 47kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human FOXP3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

JM2; AIID; IPEX; PIDX; XPID; DIETER;

Catalog #: A5706

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 50943

Isotype: IgG

Swiss Prot: Q9BZS1

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the forkhead/winged-helix family of transcriptional regulators. Defects in this gene are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX), also known as X-linked autoimmunity-immunodeficiency syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).