

GATA3

Reactivity: Human

Tested applications: WB IF FC

Recommended Dilution: WB 1:500 - 1:2000 IF 1:20 - 1:50 FC 1:20 - 1:50

Calculated MW: 48kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human GATA3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

HDR; HDRS;

Catalog #: A5711

Antibody Type:

Monoclonal Antibody

Species: Mouse

Gene ID: 2625

Isotype: IgG

Swiss Prot: P23771

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein which belongs to the GATA family of transcription factors. The protein contains two GATA-type zinc fingers and is an important regulator of T-cell development and plays an important role in endothelial cell biology. Defects in this gene are the cause of hypoparathyroidism with sensorineural deafness and renal dysplasia.

To place an order, please [Click HERE](#).