

FUT2

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 39kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human FUT2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SE; Se2; sej; SEC2; B12QTL1;

Catalog #: A5721

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2524

Isotype: IgG

Swiss Prot: Q10981

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a Golgi stack membrane protein that is involved in the creation of a precursor of the H antigen, which is required for the final step in the soluble A and B antigen synthesis pathway. This gene is one of two encoding the galactoside 2-L-fucosyltransferase enzyme. Two transcript variants encoding the same protein have been found for this gene.

To place an order, please [Click HERE](#).