

IL1R1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 65kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human IL1R1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P80; IL1R; IL1RA; CD121A; D2S1473; IL-1R-alpha;

Catalog #: A5727

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3554

Isotype: IgG

Swiss Prot: P14778

Purity: Affinity purification

For research use only.

Background:

This gene encodes a cytokine receptor that belongs to the interleukin-1 receptor family. The encoded protein is a receptor for interleukin-1 alpha, interleukin-1 beta, and interleukin-1 receptor antagonist. It is an important mediator involved in many cytokine-induced immune and inflammatory responses. This gene is located in a cluster of related cytokine receptor genes on chromosome 2q12.

To place an order, please [Click HERE](#).