

PRPF31

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 55kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PRPF31

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

RP11; PRP31; SNRNP61; NY-BR-99;

Catalog #: A5732

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26121

Isotype: IgG

Swiss Prot: Q8WWY3

Purity: Affinity purification

For research use only.

Background:

This gene encodes a component of the spliceosome complex and is one of several retinitis pigmentosa-causing genes. When the gene product is added to the spliceosome complex, activation occurs.

To place an order, please [Click HERE](#).