

ACTR2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 44kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ACTR2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ARP2;

Catalog #: A5734

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10097

Isotype: IgG

Swiss Prot: P61160

Purity: Affinity purification

For research use only.

Background:

The specific function of this gene has not yet been determined; however, the protein it encodes is known to be a major constituent of the ARP2/3 complex. This complex is located at the cell surface and is essential to cell shape and motility through lamellipodial actin assembly and protrusion. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).