

CD14

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 40kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CD14

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A5737

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 929

Isotype: IgG

Swiss Prot: P08571

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a surface antigen that is preferentially expressed on monocytes/macrophages. It cooperates with other proteins to mediate the innate immune response to bacterial lipopolysaccharide. Alternative splicing results in multiple transcript variants encoding the same protein.

To place an order, please [Click HERE](#).