

CDC7

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 63kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CDC7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Hsk1; CDC7L1; HsCDC7; huCDC7;

Catalog #: A5738

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8317

Isotype: IgG

Swiss Prot: O00311

Purity: Affinity purification

For research use only.

Background:

This gene encodes a cell division cycle protein with kinase activity that is critical for the G1/S transition. The yeast homolog is also essential for initiation of DNA replication as cell division occurs. Overexpression of this gene product may be associated with neoplastic transformation for some tumors. Multiple alternatively spliced transcript variants that encode the same protein have been detected.

To place an order, please [Click HERE](#).