

## DAPK1

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 160kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human DAPK1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

bq

**Synonym:**

DAPK;

**Catalog #:** A5741

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 1612

**Isotype:** IgG

**Swiss Prot:** P53355

**Purity:** Affinity purification

For research use only.

**Background:**

Death-associated protein kinase 1 is a positive mediator of gamma-interferon induced programmed cell death. DAPK1 encodes a structurally unique 160-kD calmodulin dependent serine-threonine kinase that carries 8 ankyrin repeats and 2 putative P-loop consensus sites. It is a tumor suppressor candidate. Alternative splicing results in multiple transcript variants.

*To place an order, please [Click HERE](#).*