

GABRR1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 55kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GABRR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

g

Catalog #: A5745

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2569

Isotype: IgG

Swiss Prot: P24046

Purity: Affinity purification

For research use only.

Background:

GABA is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA receptors, which are ligand-gated chloride channels. GABRR1 is a member of the rho subunit family. Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).