

GRIP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 122kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human GRIP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GRIP;

Catalog #: A5747

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23426

Isotype: IgG

Swiss Prot: Q9Y3R0

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the glutamate receptor interacting protein family. The encoded scaffold protein binds to and mediates the trafficking and membrane organization of a number of transmembrane proteins. Alternatively spliced transcript variants encoding different isoforms have been described.

To place an order, please [Click HERE](#).