

TCN2

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 47kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TCN2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

II; TC; TC2; TC-2; TCII; TC II; D22S676; D22S750;

Catalog #: A5755

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6948

Isotype: IgG

Swiss Prot: P20062

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the vitamin B12-binding protein family. This family of proteins, alternatively referred to as R binders, is expressed in various tissues and secretions. This plasma protein binds cobalamin and mediates the transport of cobalamin into cells. This protein and other mammalian cobalamin-binding proteins, such as transcobalamin I and gastric intrinsic factor, may have evolved by duplication of a common ancestral gene. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).