

TP53BP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 213kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TP53BP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

p202; 53BP1;

Catalog #: A5757

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7158

Isotype: IgG

Swiss Prot: Q12888

Purity: Affinity purification

For research use only.

Background:

TP53BP1 (tumor protein p53 binding protein 1) is a protein-coding gene. Go annotations related to this gene include methylated histone residue binding and p53 binding. Plays a key role in the response to DNA damage. May have a role in checkpoint signaling during mitosis. Enhances TP53-mediated transcriptional activation.

To place an order, please [Click HERE](#).