

MOAP1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 40kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human MOAP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MAP-1; PNMA4;

Catalog #: A5759

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64112

Isotype: IgG

Swiss Prot: Q96BY2

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene was identified by its interaction with apoptosis regulator BAX protein. This protein contains a Bcl-2 homology 3 (BH3)-like motif, which is required for the association with BAX. When overexpressed, this gene has been shown to mediate caspase-dependent apoptosis.

To place an order, please [Click HERE](#).