

MAVS

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:100 IP 1:20 - 1:100

Calculated MW: 56kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MAVS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

5

Synonym:

IPS1; VISA; IPS-1; CARDIF;

Catalog #: A5764

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57506

Isotype: IgG

Swiss Prot: Q7Z434

Purity: Affinity purification

For research use only.

Background:

This gene encodes an intermediary protein necessary in the virus-triggered beta interferon signaling pathways. It is required for activation of transcription factors which regulate expression of beta interferon and contributes to antiviral immunity. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).