

NR5A2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 61kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NR5A2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

B1F; CPF; FTF; B1F2; LRH1; LRH-1; FTZ-F1; hB1F-2; FTZ-F1beta;

To place an order, please [Click HERE](#).

Catalog #: A5766

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2494

Isotype: IgG

Swiss Prot: O00482

Purity: Affinity purification

For research use only.