

CDKAL1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 65kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CDKAL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A5771

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 54901

Isotype: IgG

Swiss Prot: Q5VV42

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the methylthiotransferase family. The function of this gene is not known. Genome-wide association studies have linked single nucleotide polymorphisms in an intron of this gene with susceptibility to type 2 diabetes.

To place an order, please [Click HERE](#).