

GSTO2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 28kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GSTO2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GSTO 2-2; bA127L20.1;

Catalog #: A5774

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 119391

Isotype: IgG

Swiss Prot: Q9H4Y5

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is an omega class glutathione S-transferase (GST). GSTs are involved in the metabolism of xenobiotics and carcinogens. Four transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).