

## HEXIM1

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:** 40kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human HEXIM1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CLP1; EDG1; HIS1; MAQ1;

**Catalog #:** A5775

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 10614

**Isotype:** IgG

**Swiss Prot:** O94992

**Purity:** Affinity purification

For research use only.

**Background:**

Expression of this gene is induced by hexamethylene-bis-acetamide in vascular smooth muscle cells. This gene has no introns.

*To place an order, please [Click HERE](#).*