

NDEL1

Reactivity: Human Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 38kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NDEL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

EOPA; NDE2; NUDEL; MITAP1; NDE1L1;

Catalog #: A5776

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 81565

Isotype: IgG

Swiss Prot: Q9GZM8

Purity: Affinity purification

For research use only.

Background:

This gene encodes a coiled-coil protein that plays a role in multiple processes including cytoskeletal organization, cell signaling and neuron migration, outgrowth and maintenance. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome X.

To place an order, please [Click HERE](#).