

SLAMF7

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 37kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SLAMF7

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

19A; CS1; CD319; CRACC;

To place an order, please [Click HERE](#).

Catalog #: A5782

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57823

Isotype: IgG

Swiss Prot: Q9NQ25

Purity: Affinity purification

For research use only.