

PPP2R1A

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 65kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PPP2R1A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bt

Synonym:

PR65A; PP2AAALPHA; PP2A-Aalpha;

Catalog #: A5799

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5518

Isotype: IgG

Swiss Prot: P30153

Purity: Affinity purification

For research use only.

Background:

This gene encodes a constant regulatory subunit of protein phosphatase 2. Protein phosphatase 2 is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. The constant regulatory subunit A serves as a scaffolding molecule to coordinate the assembly of the catalytic subunit and a variable regulatory B subunit. This gene encodes an alpha isoform of the constant regulatory subunit A. Alternatively spliced transcript variants have been described.

To place an order, please [Click HERE](#).