

SATB1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 85kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SATB1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A5800

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6304

Isotype: IgG

Swiss Prot: Q01826

Purity: Affinity purification

For research use only.

Background:

This gene encodes a matrix protein which binds nuclear matrix and scaffold-associating DNAs through a unique nuclear architecture. The protein recruits chromatin-remodeling factors in order to regulate chromatin structure and gene expression. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).