

SERPINA6

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 45kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SERPINA6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CBG;

Catalog #: A5803

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 866

Isotype: IgG

Swiss Prot: P08185

Purity: Affinity purification

For research use only.

Background:

This gene encodes an alpha-globulin protein with corticosteroid-binding properties. This is the major transport protein for glucorticoids and progestins in the blood of most vertebrates. The gene localizes to a chromosomal region containing several closely related serine protease inhibitors which may have evolved by duplication events.

To place an order, please [Click HERE](#).