

DEDD

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 36kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human DEDD

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DEFT; KE05; DEDD1; FLDED1; CASP8IP1;

Catalog #: A5806

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9191

Isotype: IgG

Swiss Prot: O75618

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that contains a death effector domain (DED). DED is a protein-protein interaction domain shared by adaptors, regulators and executors of the programmed cell death pathway. Overexpression of this gene was shown to induce weak apoptosis. Upon stimulation, this protein was found to translocate from cytoplasm to nucleus and colocalize with UBTF, a basal factor required for RNA polymerase I transcription, in the nucleolus. At least three transcript variants encoding the same protein have been found for this gene.

To place an order, please [Click HERE](#).