

EPC1

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 93kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human EPC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Epl1;

Catalog #: A5807

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 80314

Isotype: IgG

Swiss Prot: Q9H2F5

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the polycomb group (PcG) family. The encoded protein is a component of the NuA4 histone acetyltransferase complex and can act as both a transcriptional activator and repressor. The encoded protein has been linked to apoptosis, DNA repair, skeletal muscle differentiation, gene silencing, and adult T-cell leukemia/lymphoma. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).