

EZH1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 85kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human EZH1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

KMT6B;

Catalog #: A5818

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2145

Isotype: IgG

Swiss Prot: Q92800

Purity: Affinity purification

For research use only.

Background:

EZH1 is a component of a noncanonical Polycomb repressive complex-2 (PRC2) that mediates methylation of histone H3 (see MIM 602812) lys27 (H3K27) and functions in the maintenance of embryonic stem cell pluripotency and plasticity (Shen et al., 2008)

To place an order, please [Click HERE](#).