

HELLS

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP CHIP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200 IP 1:50 - 1:200

ChIP 1:20 - 1:100

Calculated MW: 97kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human HELLS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

LSH; PASG; SMARCA6; Nbla10143;

Catalog #: A5831

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3070

Isotype: IgG

Swiss Prot: Q9NRZ9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a lymphoid-specific helicase. Other helicases function in processes involving DNA strand separation, including replication, repair, recombination, and transcription. This protein is thought to be involved with cellular proliferation and may play a role in leukemogenesis. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).