

PRODH

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:68kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PRODH

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

POX; PIG6; HSPOX2; PRODH1; PRODH2; TP53I6;

Catalog #:A5836

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5625

Isotype:IgG

Swiss Prot:O43272

Purity:Affinity purification

For research use only.

Background:

This gene encodes a mitochondrial protein that catalyzes the first step in proline degradation.

Mutations in this gene are associated with hyperprolinemia type 1 and susceptibility to schizophrenia 4 (SCZD4). This gene is located on chromosome 22q11.21, a region which has also been associated with the contiguous gene deletion syndromes, DiGeorge and CATCH22.

Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).