

ITGB3BP

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 20kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ITGB3BP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

CENPR; NRIF3; TAP20; CENP-R; HSU37139;

Catalog #: A5839

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23421

Isotype: IgG

Swiss Prot: Q13352

Purity: Affinity purification

For research use only.

Background:

This gene encodes a transcriptional coregulator that binds to and enhances the activity of members of the nuclear receptor families, thyroid hormone receptors and retinoid X receptors. This protein also acts as a corepressor of NF-kappaB-dependent signaling. This protein induces apoptosis in breast cancer cells through a caspase 2-mediated signaling pathway. This protein is also a component of the centromere-specific histone H3 variant nucleosome associated complex (CENP-NAC) and may be involved in mitotic progression by recruiting the histone H3 variant CENP-A to the centromere. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).