

PTH

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW:12kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PTH

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

PTH1;

Catalog #:A5846

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5741

Isotype:IgG

Swiss Prot:P01270

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a hormone secreted by parathyroid cells. This hormone elevates blood Ca²⁺ level by dissolving the salts in bone and preventing their renal excretion. Defects in this gene are a cause of familial isolated hypoparathyroidism (FIH).

To place an order, please [Click HERE](#).