

RTN4R

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 50kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human RTN4R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NGR; NOGOR;

Catalog #: A5847

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 65078

Isotype: IgG

Swiss Prot: Q9BZR6

Purity: Affinity purification

For research use only.

Background:

This gene encodes the receptor for reticulon 4, oligodendrocyte myelin glycoprotein and myelin-associated glycoprotein. This receptor mediates axonal growth inhibition and may play a role in regulating axonal regeneration and plasticity in the adult central nervous system.

To place an order, please [Click HERE](#).