

SMARCAD1

Reactivity: Human Mouse

Tested applications: WB IF IP

Recommended Dilution: WB 1:500 - 1:2000 IF 1:50 - 1:100 IP 1:50 - 1:100

Calculated MW: 117kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SMARCAD1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

ETL1; HEL1; ADERM;

Catalog #: A5850

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 56916

Isotype: IgG

Swiss Prot: Q9H4L7

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the SNF subfamily of helicase proteins. The encoded protein plays a critical role in the restoration of heterochromatin organization and propagation of epigenetic patterns following DNA replication by mediating histone H3/H4 deacetylation. Mutations in this gene are associated with adermatoglyphia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

To place an order, please [Click HERE](#).