

## SRPK1

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

**Calculated MW:** 74kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human SRPK1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

b

**Synonym:**

SFRSK1;

**Catalog #:** A5854

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 6732

**Isotype:** IgG

**Swiss Prot:** Q96SB4

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a serine/arginine protein kinase specific for the SR (serine/arginine-rich domain) family of splicing factors. The protein localizes to the nucleus and the cytoplasm. It is thought to play a role in regulation of both constitutive and alternative splicing by regulating intracellular localization of splicing factors. Alternative splicing of this gene results in multiple transcript variants. Additional alternatively spliced transcript variants have been described for this gene, but their full length nature have not been determined.

*To place an order, please [Click HERE](#).*