

NPRL2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 43kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NPRL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

NPR2; NPR2L; TUSC4;

To place an order, please [Click HERE](#).

Catalog #: A5859

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10641

Isotype: IgG

Swiss Prot: Q8WTW4

Purity: Affinity purification

For research use only.