

TWF2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 39kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TWF2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

A6r; A6RP; PTK9L; MSTP011;

Catalog #: A5860

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11344

Isotype: IgG

Swiss Prot: Q6IBS0

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene was identified by its interaction with the catalytic domain of protein kinase C-zeta. The encoded protein contains an actin-binding site and an ATP-binding site. It is most closely related to twinfilin (PTK9), a conserved actin monomer-binding protein.

To place an order, please [Click HERE](#).