

ZBTB16

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 74kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ZBTB16

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

PLZF; ZNF145;

Catalog #: A5863

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7704

Isotype: IgG

Swiss Prot: Q05516

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the Krueppel C2H2-type zinc-finger protein family and encodes a zinc finger transcription factor that contains nine Kruppel-type zinc finger domains at the carboxyl terminus. This protein is located in the nucleus, is involved in cell cycle progression, and interacts with a histone deacetylase. Specific instances of aberrant gene rearrangement at this locus have been associated with acute promyelocytic leukemia (APL). Alternate transcriptional splice variants have been characterized.

To place an order, please [Click HERE](#).