

TFRC

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 84kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TFRC

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

T9; TR; TFR; p90; CD71; TFR1; TRFR;

To place an order, please [Click HERE](#).

Catalog #: A5865

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7037

Isotype: IgG

Swiss Prot: P02786

Purity: Affinity purification

For research use only.